


E200ix **Trauma** Stretcher




# Specifications

## Cogent E200ix Trauma Stretcher

*This stretcher with entire X-Ray function designed for whole body. It also added 90 degrees sitting posture and another function is that installed the pneumatic cylinder to adjust crook for knees from 0 degree to 40 degree.*

- Integrated storage tray with oxygen bottle trough
- Trend / Reverse Trend
- Control hydraulic pump in dual-side foot pedal
- X Ray translucent
- Restraint straps with velcro
- Brake pedal control four wheels / steering pedal in two end og bed
- IV pole inserted hole x 5
- Protective corner bumpers x4
- Collapsible siderails
- Pneumatic assist backrest bends the angle  $0^{\circ} \sim 90^{\circ}$
- Retractable fifth-wheel steering
- 4"(10cm) mattress or more
- (tente) 8" wide wheel
- Wheel Cover
- Delivery Push Bar
- Two stages movable IV pole x 1
- Two stages permanent folding IV Rod x1
- Foot auxiliary push bar
- No-crank pneumatic knee gatch  $0^{\circ} \sim 40^{\circ}$
- Upright oxygen holder aperture



# Specifications

Model Number	ES280 / E200 / E200i / E220LED A.B	E200X / E200ix	ES280
Overall Length	208 cm	212cm	Electric Specification
Overall Width	85 cm	85cm	Voltage AC220-240 / 120VAC 50/60 Hz
Weight Capacity	225 kg	228kg	Standard spare parts
Height Range			Mattress, storage battery and charger
High	91 cm ±1	99 cm ±1	
Low	55 cm ±1	65 cm ±1	

Litter Positioning	ES280 / E200 / E200i / E220LED B	E200X / E200ix
Backrest	0-90° ±1°	
Trend. / Reverse Trend	±18°	
Patient Surface	195 x 66 cm	195 x 66 cm
Footboard lift crank	0-40° ±2° ( ES280 / E200i / E220 LED. A )	
Siderails	126 x 42 cm	126 x 39 cm
Caster Diameter	8" ( 20 cm ) TENTE	8" ( 20 cm ) TENTE
No-crank pneumatic knee gatch	0-40° ( E200ix )	

## Optional Accessory

- ✓ Upright oxygen bottle trough
- ✓ Retractable fifth-wheel Steering
- ✓ Upright chest X-RAY frame
- ✓ 4"10 cm mattress or more
- ✓ Patient buckle up restraint straps
- ✓ Patient restraint straps with velcro
- ✓ Instrument holder with monitor or ventilator
- ✓ Decorative strips of three optional colors of red, sky blue and navy blue.
- ✓ Four-Side brake & steering pedals
- ✓ LED.Caution light
- ✓ Head / Foot board for bed