

Haier







-86°C ULT Freezer



DW-86L490

ULT Freezer

This product is specially designed and manufactured for long term storage of various biological products, including viruses, germs, erythrocytes, leucocytes and cutis. Applications can be found in research institutes, laboratories in electronic and chemical plants, biological engineering institutes and marine fishery companies.

	<p>Energy-efficient Double Door</p> <p>Unique double door to store frequently used samples in the top chamber while long-term storage remains undisturbed in bottom</p>		<p>Excellent Door Seals</p> <p>provides maximum protection against thermal intrusion</p>		<p>Control panel</p> <p>Microprocessor control,digital display, adjustable temperature range :-20°C--40°C.</p>
	<p>Individual Inner Doors</p> <p>4 individual inner doors can be opened independently to minimize frost buildup inside the chamber</p>		<p>Safety</p> <p>Multiple malfunction alarms include high/low temperature, power failure, sensor error, low battery, high/low ambient temperature.</p>		<p>Optional chart recorder</p> <p>Built-in chart recorder can be installed at eye level</p>

Electrical Data	Power Supply(V/HZ)	220-240/50	220/60
	Power(W)	1100	1100
	Electrical Current(A)	6.5	6.5
Dimensions	Capacity(L/Cu.Ft)	490/17.3	
	Net/Gross Weight (approx)	295/335(kg) 650.4/738.5(lbs)	
	Interior Dimensions (W*D*H)	590*630*1310(mm) 23.2*24.8*51.6(in)	
	Exterior Dimensions (W*D*H)	860*900*1980(mm) 33.9*35.4*78.0(in)	
	Packing Dimensions (W*D*H)	920*990*2165(mm) 36.2*38.9*85.2(in)	
	Container load (20'/40'/40'H)	12/24/24	
Accessories	Caster	Y	
	Foot	Y	
	Test Hole	Y	
	Shelves/Inner doors	3/4	
	Temp Recorder	Optional	

Haier Medical and Laboratory Products Co., Ltd

Room 403D, Brand Building, Haier industrial park, No. 1 Haier Road, Qingdao,266101,P.R.China
Tel: +86-0532-88936011 / 7169 Website: www.haiermedical.com

